

# How to Collect in ProteX



Consult with your fertility team for specific media and lubricant usage instructions.

- 1. Open** your sterile sealed package and place ProteX™ on a flat, stable surface.

If using refrigerated media, ensure it comes to room temperature first.

- 2. Pour** all media from vial carefully into the ProteX.

- 3. Collect** your semen sample into the ProteX container as directed by your fertility team.

If using lubrication, be sure it has been cleared by your fertility team and be careful not to get the lubricant in the ProteX.

- 4. Secure** the lid carefully onto the ProteX container.

- 5. Return** your sample to your fertility team as directed; keep ProteX upright while transporting.